


# Urogenx™ Bovine Pericardial Tissue Patch


## Urology Indications

- Cervical anastomotic leakage
- Corporoplasty
- Partial nephrectomy
- Pelvic Organ Prolapse Repair
- Peyronie's repair
- Rectovaginal fistula
- Urethroplasty
- Vesicovaginal Fistula Repair



- Artificial Urinary Diversion
- Artificial Urinary Sphincter Placement
- Bladder Augmentation
- Hypospadias Repair
- Pubovaginal Sling Surgery
- Ureteral Reconstruction
- Urethral stricture Repair
- Urologic Venous Reconstruction

## elixP™ fixated patches excel in the 4 major categories



Ultimate Tensile Strength



Suture Retention



Burst Strength



Elasticity & Elongation

## Wide Variety of Sizes with Uniformity

<b>Model</b>	<b>XM-05*</b>	<b>XM-07*</b>	<b>XM-08*</b>	<b>XM-12</b>	<b>XM-14</b>	<b>XM-15</b>
<b>Size (cm)</b>	0.8 x 8	1 x 10	1 x 14	2 x 9	4 x 4	4 x 6
<b>Model</b>	<b>XM-18</b>	<b>XM-19</b>	<b>XM-20</b>	<b>XM-21</b>	<b>XM-22</b>	<b>XM-23</b>
<b>Size (cm)</b>	5 x 10	6 x 8	6 x 10	7 x 10	8 x 14	10 x 16

\*Also available in tapered size.



Call the number or visit our website below to place an order for our products or speak to our sales representatives for more information!

+1 (949) 670-0403  
info@tisgenx.com  
www.tisgenx.com



Tisgenx, Inc.  
15615 Alton Parkway,  
Irvine, CA 92618, USA

- Bhende, Vishal V et al. "The Utility of Invengenx® Bovine Patch for Right Ventricular Outflow Tract (RVOT) Reconstruction and Augmentation in the Surgical Management of Tetralogy of Fallot (TOF): A Contemporary Study and Review of the Literature." *Cureus* vol. 15,10 e46882. 12 Oct. 2023, doi:10.7759/cureus.46882
- Bhende, Vishal V et al. "Beyond Synthetics: Promising Outcomes With the Invengenx® Bovine Pericardial Patch for Ventricular Septal Defect Repair in a Young Pediatric Population." *Cureus* vol. 16,3 e55530. 4 Mar. 2024, doi:10.7759/cureus.55530
- Shah, P., T. Dudde, S. Kamath, K. Surya, K. V. Bilgi, and R. Kupumbati. "Initial Experience of Using the Invengenx® Bovine Pericardial Tissue Patch for Common Arterial Trunk Repair". *International Journal of Research in Medical Sciences*, vol. 12, no. 7, June 2024, pp. 2722-3, doi:10.18203/2320-6012.ijrms20241941.
- Shah P, Ramasami H, Kamath S, Surya K, Bilgi KV, et al. (2024) Double Patch Repair for Complete Atrioventricular Septal Defect Using the Invengenx® Patch. *Int J Cardiol Res* 13:2.
- Weiss, S., et al. "Self Made Xenopericardial Aortic Tubes to Treat Native and Aortic Graft Infections." *Journal of Vascular Surgery*, vol. 66, no. 6, 2017, p. 1914, doi:10.1016/j.jvs.2017.10.007.
- Wiggins, Luke M., et al. "The Utility of Aortic Valve Leaflet Reconstruction Techniques in Children and Young Adults." *The Journal of Thoracic and Cardiovascular Surgery*, vol. 159, no. 6, 2020, pp. 2369-2378, doi:10.1016/j.jtcvs.2019.09.176.


## elixP™ Tissue

### What is elixP™?

At the heart of our innovative bovine pericardial tissue patch lies elixP™, a cutting-edge processing technique that sets a new standard in medical technology. Our proprietary elixP™ Fixation Technology is meticulously engineered to preserve the natural triple helical structure of collagen at both intramolecular and intermolecular levels. This breakthrough process ensures that the collagen molecules retain their intrinsic alignment and integrity, resulting in a tissue patch with superior mechanical strength and flexibility.


### elixP™ tissue technology

elixP™ Tissue Technology achieves a remarkable level of precision, with 100% crosslinking of collagen molecules. This unparalleled level of crosslinking is critical for ensuring the patch's stability and longevity. By solidifying the bonds within and between collagen fibrils, elixP™ effectively prevents issues such as suture line bleeding, delamination, and inflammatory responses. The result is a tissue patch that remains intact and functional throughout the healing process, minimizing potential complications and improving surgical outcomes. With elixP™, we eliminate antigenicity and maintain the natural collagen formation, providing a patch that boasts exceptional softness and durability.




### Clinical benefits and long-term results

Our advanced elixP™ tissue technology not only supports superior biocompatibility but also resists degradation, allowing Urogenx™ to perform reliably across a range of surgical applications. The elixP™ methodology does more than just crosslink collagen—it transforms the collagen matrix into a superior material that is both highly durable and exceptionally flexible. This enhanced matrix allows for precise manipulation during application, offering excellent scaffolding and maintaining the natural pliability of the tissue. Surgeons can rely on Urogenx™ for its ease of application and adaptability, knowing that it will conform perfectly to the surgical site and support effective tissue repair. The enhanced tensile strength and suture-holding capabilities of Urogenx™ further contribute to its reliability and effectiveness. Urogenx™ represents a significant advancement in tissue repair technology, offering a solution that medical professionals can trust for complex and critical procedures.



**elixP™ tissue under H&E microscopy.**  
Collagen pattern is almost identical to native pericardium.



**Competitor tissue under H&E microscopy.**  
Large white gaps between collagen leads to infection and calcification.